



**CLARITY MEDICAL**

claritymedical.healthcare

Clarity Medical Ltd  
t. +44 (0)330 124 7287

## G1 Datasheet

# Natural Latex Powder Free Gloves (G1)

### General Details

Material	Natural Latex
Surface Finish	Powder Free (inner chlorinated, outer polymer coated)
Colour	Natural
Geometry	Ambidextrous
Texture	Fully Textured (Palm Textured)
Cuff and Finishing	Beaded

### Glove Specifications

Size	Weight g/pc	Length mm	Palm Width mm	Thickness	
				Cuff	Palm
XS	4.8 ± 0.2	240	≤ 80	0.09 ± 0.01	0.13 ± 0.01
S	5.3 ± 0.2	240	85 ± 5	0.09 ± 0.01	0.13 ± 0.01
M	5.8 ± 0.2	240	95 ± 5	0.09 ± 0.01	0.13 ± 0.01
L	6.4 ± 0.2	240	105 ± 5	0.09 ± 0.01	0.13 ± 0.01
XL	6.9 ± 0.2	240	≥ 110	0.09 ± 0.01	0.13 ± 0.01
Powder Content (Maximum)			2 mg / glove		
Powder Content (Maximum)			50 µg / g		

AQL 1.5  
EN455 Parts 1, 2, 3 & 4  
MDD 93/42/EEC  
EN374

Class 1 Medical Device



### Performance Requirements

Test	Related Defects	Inspection Level	AQL
Water Leak at Final Release	Holes	G-1 Single Sampling Plan	1.5

### Manufacturing Standards

Manufacturing Standards	Incompliance with ISO 9001:2008
Presentation	Standard: 100 gloves per dispenser box; 10 dispenser boxes per Carton
Glove Care	Should not store under direct sunlight or at high heat and humid conditions. Single use only.